


尺寸:
70*100mm



eSebone

Multi-action Full Back & Car Seat massage

USER MANUAL



HC 1048

Design & Developed in Japan

PRODUCT FEATURES


This intelligent massage cushion features 6 massage nodes to provide kneading massage on your shoulder, back and waist; Vibration massage for buttocks with 2 intensity settings; Back massage nodes travel up and down with 3 modes: full back, upper or lower, you may also focus massage on specific back areas according to your liking; Optional soothing heat for back helps relieve fatigue and loosen muscle stiffness; 15-minute auto-off function and overheat protection design promote safer use; Easy-to-use controller allows quick control of all functions.

OPERATING INSTRUCTION


Preparation before use

-Check whether the product is damaged or not
-place the product on the body part to be massaged, and make sure that the place is appropriate.

Methods Connect the DC head of the power adapter into the DC line interface of the product, and insert the other end into the power socket to start the massage through the hand controller. It defaults to a 15 minutes massage session. Press this button to switch to 25 minutes, press again to switch to 10 minutes.



BUTTON DESCRIPTION



1 Power On/Off
Press this button to start, back massage nodes travel up and down delivering kneading massage along the full back, you may press kneading button to switch to upper back or lower back, or press point button to focus the massage on current point. Back heat is activated, seat cushion also starts low speed vibration.

Heat
When massage cushion is working, press this button to turn on/off back heat function.

Vibration
When massage cushion is working, press this button to switch to high speed vibration, press again to turn off seat vibration function.

Kneading
It defaults to a full back kneading massage. When full back kneading massage is running, press this button to switch to upper back, press again to switch to lower back, press for the third time to stop back kneading massage.

Point
When back kneading massage is working, press this button to focus massage to current point, press again to stop point massage and switch to full back kneading massage.

Time
It defaults to a 15 minutes massage session. Press this button to switch to 25 minutes, press again to switch to 10 minutes.



SPECIFICATION DATA OF PRODUCTS

Type: Hc1048
Power: 24W
Rate voltage: 100-220V~
Rate frequency: 12V ---
Working time: 10/15/25 minutes

APPLICABLE PEOPLE

1. People who use mobile phones for a long time
2. People who work with computers for a long time.
3. People who are nervous, insomnia and dreamy.
4. People with stiff shoulders and neck and muscle fatigue
5. People with sore shoulders and neck and weak hands.

TROUBLE CLEARING

Symptoms	Possible Causes	Solution
Appliance doesn't start	Power strips is not working	Check the power strip to see if it is good Check the connection between power plug and power strip
Power strip is good but appliance still not working	Appliance is faulty	Unplug the appliance from the electrical outlet and contact front service team



STORAGE AND MAINTENANCE

1. Can't be close to the ground, walls and roofs.
2. The room should be well ventilated to avoid strong sunlight and corrosive odor.
3. The product is equipped with simple anti-vibration facilities, which is suitable for aviation, railway, highway and ship transportation. Rain and snow splash, inversion and collision should be avoided.

Maintenance

1. It must be moistureproof, rat proof and pollution proof.
2. When cleaning the handheld controller, it should be done in the state of shutdown.
3. After use, the product should be placed in a safe place. Do not put heavy objects on the product.
4. When something goes wrong, please ask professional technicians to repair it, and don't disassemble the product at will.

MATTERS NEED ATTENTION

1. Do not use massager while driving.
2. Do not use the massager in high temperature and humidity.
3. Before using the massager, please bundle up your hair for those with long hair.
4. For those suffering from hypertension and mental illness, please do not use this product when they are ill.
5. This massager is a personal care product, which can't replace the advice and treatment provided by medical professionals. If you have any questions or needs, please consult your medical professionals.
6. This massager is designed as a household electronic health care product, which is not suitable for diagnosis, treatment and other commercial purposes.
7. If cleaning is required, it should be wiped with sterilized cotton soaked in neutral detergent. Do not wash the product with water.



WARRANTY METHOD

The product is guaranteed free of charge within one year from the date of purchase. Please show the warranty card when repairing. If the product has functional problems since the date of sale, you can return/exchange it at the purchase counter with the product warranty card or invoice within 7 days. The return/exchange is not applicable in any circumstances. If the product has functional problems within 15 days from the date of sale, you can choose to exchange or repair it.

COMPANY INFORMATION



1. Your HCL product is fully under warranty against servicing of any manufacturing defects for the period of 12 months from the date of purchase unless stated otherwise.
2. Manufacturing defect towards any of the spare parts of any machinery will be covered for Six Months from the date of Purchase. Upholstery is not covered under the warranty.
3. This warranty is not transferable neither exchange with any other product. It is only applicable to original purchaser and in the city of purchase.
4. This warranty shall be subject to the following limitations and exclusions: -
(a) Malfunction or damages resulting from natural disaster, fire, civil unrest and or accidents.
(b) Defects from using wrong electrical supply, Misuses or damages by negligence are and abuse.
() Accessories.
5. With HCL wide network you send purchase for servicing to any HCL service centre around the world. If that is not the same place of purchase, Individual Country, State, or City services policy applies.
6. To find the HCL centre please contact the HCL representative. In case of any fault the product must be brought to the service centre. The customer will bear the cost of transportation to and from the service centre.
7. Sales made to cities/Places outside the operations city will also be entitled to have warranty from the nearest operational city till such time service centre is opened in their city. However, customer will bear the cost of transportation to and from Service centre.
8. It is advisable to use CVT to avoid voltage fluctuations.

Sales Office:
Service Support: +91-73000-00739

Customer Name:	Invoice No:
Address:	
Pin:	State:
Phone Number:	Cell phone number:
Email ID:	Date:
Country:	

Specifications:
Product: eSebone
Model No. HC 1048
Voltage: 100-240W (DC 12V, 2A)
Power: 30W
Auto Time: 15 M
D&D: Japan

Sales & Service Inquiry:
+91-73000-00739
www.humanicare.com
info@humanicare.com